

Arrows to indicate the opening direction

For expansion screws.

Color	REF	Quantity
yellow	611-130-00	50 pieces
white	611-131-00	50 pieces
pearl-blue	611-132-00	50 pieces



CE

Universal bending key

Easy bending of the retention legs of all hyrax® and Variety screws without using any additional pliers (leverage principle).
Material: stainless steel

REF 611-122-00	1 piece
----------------	---------

**Safety key**

Can be used with all Dentaaurum hyrax® and Variety screws.
For safe intraoral activation of the expansion screws in fixed appliances.

REF 611-120-01	10 pieces
----------------	-----------

CE

**Expansion screw key
straight, small – standard**

Can be used with all Dentaaurum expansion screws.
To activate all expansion screws in removable appliances (Outside of the mouth).
Material: stainless steel

REF 611-115-00	1 piece
REF 611-115-30	10 pieces
REF 611-115-60	100 pieces

**Expansion screw key
large – with elongated handle**

Can be used with all Dentaaurum expansion screws.
To activate all expansion screws in removable appliances (Outside of the mouth).
Material: stainless steel

REF 611-116-00	1 piece
----------------	---------

**Expansion screw sample kit**

Filled with 32 expansion screws that are also available separately, 2 expansion screw keys and 1 universal bending key. **1 set = 35 pieces.**

REF 612-000-00	1 set
----------------	-------

CE 0483

