

Extraoral appliances



High-pull headgear for chin cap therapy

Suitable for all head and chin sizes.

The connector line runs between the tip of the chin and temporomandibular joint.

Indications:

- for early treatment of Angle Class III
- with cross-bite tendency of single or all anterior teeth

Chin cap	Color	REF	Quantity
elastic, porous and air permeable	blue	744-705-00	1 piece
rigid, short outer hooks, air permeable	blue	744-709-00	1 piece

CE



Vertical-pull headgear for chin cap therapy

Suitable for all head and chin sizes.

The connector line runs between the tip of the chin and temporal region.

Indications:

- for early treatment of an open bite
- for early treatment of a vertical growth pattern

Chin cap	Color	REF	Quantity
elastic, porous and air permeable	blue	744-706-00	1 piece
rigid, short outer hooks, air permeable	blue	744-710-00	1 piece

CE



High-pull headgear "Hickham"

Suitable for all head sizes.

The required force is produced by extraoral elastics (see accessories).

Indications:

- for treatment of Angle Class II/2 (anterior deep overbite)
- for intrusion of anterior teeth
- for distalization of cuspids/bicuspid

J-Hooks (REF 744-713-00) for Hickham headgear see accessories.

Color	REF	Quantity
blue	744-717-00	1 piece

CE