


## Instruments



**Molar band seater**  
**sterilizable to a maximum of 140 °C / 284 °F**  
 With replaceable plastic head and hard metal triangular tip.

REF 026-354-00	1 piece
----------------	---------

CE



**Plastic head with triangular hard metal tip**  
**sterilizable to a maximum of 140 °C / 284 °F**  
 For molar band seater, REF 026-354-00.

REF 026-357-00	1 piece
----------------	---------

CE



**Mershon band pusher**  
 acc. to Mershon  
**sterilizable**  
 With serrated rectangular tip and anatomically formed grip for good transfer of force when positioning bands.

REF 024-300-00	1 piece
----------------	---------

CE



**Bracket positioning gauge**  
**sterilizable**  
 Simplifies seating the brackets on the teeth in the correct distance between bracket slot and incisal edge.  
**With metal tips.**

REF 030-390-00	1 piece
----------------	---------

CE



**Bracket positioning gauge**  
 Simplifies seating the brackets on the teeth in the correct distance between bracket slot and incisal edge.  
**With pencils.**  
 Not approved for intraoral use!

REF 030-391-00	1 piece
----------------	---------



**Safety caps**  
 For bracket positioning gauge, REF 030-391-00.  
 Not approved for intraoral use!  
 Color: white

REF 030-392-00	4 pieces
----------------	----------